



Lucira ankaptik COVID-19 Test Kit pitek ta'n keknuwa'tikek

Kite'n msit ta'n ewikasik ke'smna'q poqtamkiaq


<p>Test kit ta'n koqowey pitek</p> 	<p>Ke'smna'q poqtamkianuk</p>  <p>Waqma'ten kpitnn</p> <p>Wi'ke'n wisun test- ituk</p>
---	---

Lukwatis ta'n saqte'k. Kwile'n ta'n ma wen like aqq atestmulnin ta'n tujiw eskman 30 minutes test. Tlwula'tike aqq pipanikesi ta'n tujiw kitnman majukwatmn ta'n teluwek. Jiptuk ki'l tlo'tis l'nu'jaqmati wela'sik.

1 Ikanluwate'n wjit ilajuktmin test

 <p>Panta'te'n packages 1 aqq 2</p>	 <p>Pija'te'n tapu'kl AA batteries</p>	<p>Amuj pa saqte'k</p>  <p>Ankapte'n stoqnamu'kn wasoqeqkn</p>	 <p>Mna'te'n seal</p>	<p>menaqaj</p>  <p>Mukk paskija'tu</p>	 <p>Mukk panta'tu swap misoqo kuwayjk e'wmn</p>
--	--	--	--	--	--

2 Sisqon ika'te'n swab

 <p>Mukk sama'tu</p> <p>Mukk paskija'tu</p>	 <p>1 2</p> <p>2cm 2cm</p> <p>5x 5x</p> <p>Kasipute'n te's sisqon etu'k</p> <p>Msit swab pija'tu'n misoqo mist pitek ksisqonk Jiptuk tl'otis wikawe'k.</p>
---	--

3 Kito'qa'te'n swab vial

 <p>Siawa'te'n ta'n metoqiaq</p>	 <p>15x</p>	<p>Jikleketisk swab</p> 	 <p>Siawa'sik piluwey page</p>
---	--	--	---

4 Ankapte'n wula'sin test



KEKNUWA'TASIK wasoqitesk te's 5 seconds. Katu mu wesoqatesknuk, paskija'te'n ma'minu vial.

Skmate'n 30 minutes



Pasik MU aja'tu test



Ta'n tujiw kaqkisa'tu'n, la'tu ta'n maliaptik airline wjit ta'n teluwek



Majukwate'n airline ta'n elukutijik telimiskik wjit ta'n teli jikleken test kit

Keknuwa'tikek

? Mu kelu'knuk kisa'tasik



Pipanikesites piluwey kit aqq apa suwa'tis test

- Teluwek mu weskunumu'n



Nike' kisi siawa'sitis tujiw lietis ta'n pemilien Ula teluwek pasik eliklu'k 72 hours

+ Teluwek weskunmin



COVID-19 we'jitasik. Mukk we'kwata'siw, aqq knuwa'tu ta'n ikana'toq airline amkmayu. Jiptuk ma kisi liwun ta'n pemi lien.

Jiptuk amuj pa majukwatis wejiaq airport staff, airline staff, aqq ta'n pema'tu'tij wlo'ti wjit ta'n tlukwetis.